




Alexis[®]

Orthopaedic Protector



Alexis


Orthopaedic Protector

Shields soft tissue – including skin, fat, and muscles – from contact with abrasive instrumentation and debris at the wound margins.


Distributes the force of retraction evenly for excellent exposure of the joint capsule.

Provides 360 degrees of atraumatic protection and retraction.

Total Hip Arthroplasty



Total Shoulder Arthroplasty





Alexis Orthopaedic Protector

Featuring a rigid retraction ring for maximum retraction and a flexible retraction ring for maximum versatility*

Reorder No.	Size	Sheath Length	Incision Range	Qty/Box
-------------	------	---------------	----------------	---------

Rigid Retraction Ring

HR000	X-Small/Medium	14cm	2.5–7cm	5
HR001	Small/Small	14cm	2.5–8cm	5
HR004	Small/Medium	14cm	2.5–8cm	5
HR005	Medium/Large	17cm	5–13cm	5

Flexible Retraction Ring

HR100	X-Small/Medium	14cm	2.5–7cm	5
HR101	Small/Small	14cm	2.5–8cm	5
HR104	Small/Medium	14cm	2.5–8cm	5
HR105	Medium/Large	17cm	5–13cm	5

*Compared to metal retractors.

Visit [appliedmedical.com/alexis](https://www.appliedmedical.com/alexis) for more information.

The Alexis orthopaedic protector is used to access the musculoskeletal system through an atraumatically retracted wound and shield soft tissue from contact with abrasive instrumentation and debris.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

09/2022 431289-EN-C

