



Inzii[®]

Retrieval Systems

A large teal graphic element consisting of two overlapping diagonal shapes that form a triangular shape pointing towards the bottom right corner of the page.

COMPLETE SPECIMEN RETRIEVAL BAG OFFERING

Reinforced tip design ensures specimen security and
reduces risk of contamination

Innovative design allows introduction through smaller
ports, advancing minimally invasive procedures


Bag remains open and intact throughout
specimen capture



Inzii Universal 10mm bag in a 5mm* device

Reduces the need for large ports and enhances versatility

Extra-large angled design eases specimen delivery into bag





Inzii 12/15mm

15mm bag in a 12mm* device



Inzii 10mm

Standard specimen retrieval bag

Procedural Applications


Facilitates a wide range of laparoscopic procedures, including:



LAPAROSCOPIC

Cholecystectomy
Oophorectomy
Appendectomy

Nephrectomy
Splenoectomy
Sleeve Gastrectomy



Based on an unpublished internal test report from September 2013 and November 2015

Guide Bead

Allows reopening and closure of bag, facilitating multiple specimen retrieval





Inzii Universal Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD003	5mm*	180ml	10/Box



Inzii 10mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD001	10mm	225ml	10/Box



Inzii 12/15mm Retrieval System

Reorder No.	Port Size	Volume	Qty/Box
CD004	12mm*	1600ml	5/Box

* Applied Medical Kii Access System

Visit [appliedmedical.com/inzii](https://www.appliedmedical.com/inzii) for more information.

The Inzii retrieval systems are used during procedures for specimen retrieval.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

06/2022 451624-EN-A

