



Inzii™

Retrieval Systems

A VERSATILE OFFERING OF SPECIMEN RETRIEVAL BAGS

Inzii retrieval systems, with bags available in either polyurethane material or ripstop nylon material, facilitate specimen extraction through small linear defects. Offering clinical versatility, durability and ease of use, Inzii retrieval systems maximize specimen security during tissue removal.



**Inzii Retrieval Systems –
Polyurethane Bags In-Service**



**Inzii Retrieval Systems –
Ripstop Nylon Bags In-Service**

Inzii


Retrieval Systems (Polyurethane Bags)

All Inzii polyurethane bags are 4 mils thick and have a 6-mil reinforced tip, ensuring specimen security and minimizing risk of contamination.



Inzii universal retrieval system

Delivers a 10mm bag through a 5mm port,* enhancing clinical versatility and reducing the need for large ports.



Inzii 10mm retrieval system

Provides unmatched strength and durability, as shown in the chart below.




Inzii 12/15mm retrieval system

Delivers a 15mm bag through a 12mm port, enhancing clinical versatility and reducing the need for large ports.

* Compatible with Applied Medical Kii™ 5mm trocars and larger trocars.

Burst Strength Comparison



Note: Data is on file.

Inzii

Retrieval Systems (Ripstop Nylon Bags)

Inzii ripstop nylon bags provide maximum specimen security and are ideal for cases that require increased bag strength.



Inzii ripstop redeployable retrieval system

Allows for multiple deployments and multiple specimen extractions within the same procedure, using ripstop nylon for bag durability.



Inzii 12/15mm ripstop retrieval system

Delivers a durable, 15mm ripstop bag through a 12mm port, enhancing clinical versatility and reducing the need for large ports.

Inzii

Retrieval Systems – Features




Strong, Durable Bags

Ensure secure specimen extraction.



Unique Guide Bead

Makes it easy to reopen and close the bag.



Wide Opening

Allows for smooth delivery of the specimen into the bag.

Polyurethane Bags

Description	Reorder No.	Port Size	Volume	Qty/Box
Inzii universal retrieval system	CD003	5mm ^a	180mL	10
Inzii 10mm retrieval system	CD001	10mm	225mL	10
Inzii 12/15mm retrieval system	CD004	12mm	1,600mL	5

^a Model CD003 is compatible with Applied Medical Kii 5mm trocars and larger trocars.

Ripstop Nylon Bags

Description	Reorder No.	Port Size	Volume	Qty/box
Inzii ripstop redeployable retrieval system	CD005	11mm	265mL	5
Inzii 12/15mm ripstop retrieval system	CD008	12mm	1,600mL	5

Visit [appliedmedical.com/inzii](https://www.appliedmedical.com/inzii) for more information.

Inzii retrieval systems are used during procedures for specimen retrieval.

Devices listed must be administered by a healthcare professional. Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2025 Applied Medical Resources Corporation, Rancho Santa Margarita, California. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol TM are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

All other product names, logos, brands, and registered or unregistered trademarks referenced herein are the property of their respective owners, and no trademark rights to the same are claimed. Use of these names, trademarks and brands does not imply endorsement.

