



# Applied<sup>®</sup> Scope Warmer



# A SIMPLE SOLUTION FOR BETTER VISIBILITY

## The Answer to Laparoscopic Fogging

The Applied scope warmer retains heat throughout lengthy procedures, eliminates the need for chemical gels or wipes, and minimizes the need for repeated scope cleansing and defogging.

### Scope-Warmer Seal

Retains heat throughout lengthy procedures.

### Scope-Warmer Base

Stabilizes the system when a scope is inserted.

### Scope-Warmer Bottle

Minimizes fogging by prewarming the laparoscope to body temperature.





## Applied Scope Warmer

Reorder No.	Description	Qty/Box
C3001	Scope-warmer bottle* Reusable Volume capacity = 500mL	1
C3002	Scope-warmer base* Reusable	1
C3101	Scope-warmer seals Disposable, single use Individually packaged in sterile pouches	20

\* Recommended sterilization method: Autoclave or STERRAD\*\*100S system. Please see the Instructions for Use (IFU) for complete information. The scope is not included.

Visit [appliedmedical.com/laparoscopic](https://www.appliedmedical.com/laparoscopic) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2021 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol © are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

11/2021 411033-EN-USA-A

