




Applied[®] Insufflation Needle

ENABLING HIGH-FLOW INSUFFLATION

13-gauge percutaneous needle
Tactile and visual indicators
Ergonomic housing design

Applied Insufflation Needle



Spring-Loaded, Blunt Stylet with Red Indicator Button

Provides tactile and visual feedback when the stylet is retracted.

13-Gauge High-Flow Needle

Offers faster insufflation than traditional needles, reducing operating room time.



Ergonomic Housing Design

Ensures control during needle insertion.





Applied Insufflation Needle

Reorder No.	Size	Qty/Box
C2201	13 gauge, 120mm	20
C2202	13 gauge, 150mm	20

Visit [appliedmedical.com/laparoscopic](https://www.appliedmedical.com/laparoscopic) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

03/2022 411136-EN-USA-B

