




AeroStat[®]

Cholangiography Catheter

ONE HAND ONE STEP

Easy duct access
Self-anchoring
One-handed activation



Easy Duct Access

Retention cone eliminates problems traditionally associated with laparoscopic cholangiography, such as clipping, balloon purging and potential breakage.

Self-Anchoring

Soft 6F-20F expandable cone and 3F tapered rigid tip atraumatically engages and occludes the cystic duct, eliminating the need for clips or graspers.

One-Handed Activation

Ratcheting handle, with one-handed activation, makes placement fast and easy.





AeroStat Cholangiography Catheter

Reorder No.	Description	Shaft Diameter	Length
C1002	Rigid Tip with Percutaneous Introducer	6F	33cm

Visit [appliedmedical.com/laparoscopic](https://www.appliedmedical.com/laparoscopic) for more information.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2021 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol © are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

