




# AeroStat®

Cholangiography Catheter

# ONE HAND ONE STEP

Easy duct access  
Self-anchoring  
One-handed activation



## Easy Duct Access

Retention cone eliminates problems traditionally associated with laparoscopic cholangiography, such as clipping, balloon purging and potential breakage.

## Self-Anchoring

Soft 6F-20F expandable cone and 3F tapered rigid tip atraumatically engages and occludes the cystic duct, eliminating the need for clips or graspers.

## One-Handed Activation

Ratcheting handle, with one-handed activation, makes placement fast and easy.





## AeroStat Cholangiography Catheter

Reorder No.	Description	Shaft Diameter	Length
C1002	Rigid Tip with Percutaneous Introducer	6F	33cm

Visit [www.appliedmedical.com/laparoscopic](http://www.appliedmedical.com/laparoscopic)  
for more information

© 2018 Applied Medical Resources Corporation. All rights reserved.  
Applied Medical, the Applied Medical logo design and marks designated with a ® are  
trademarks of Applied Medical Resources Corporation, registered in one or more of the following  
countries: Australia, Canada, Japan, South Korea, the United States, and/or the European Union.  
SC01096C-EN

