

APPLIED MEDICAL PRODUCT CATALOG

 APPLIED MEDICAL PRODUCT CATALOG



For ordering information:
[AppliedMedical.com/ContactUs](https://www.appliedmedical.com/contact-us)
196381-EN



APPLIED MEDICAL is dedicated to developing and providing innovative products that improve patient outcomes and enable the advancement of minimally invasive surgery.

LATEX STATEMENT

Due to our unwavering commitment to patient care, Applied Medical products do not contain natural rubber latex.



APPLIED MEDICAL TRADEMARKS

The fact that a name or logo does not appear on this list does not constitute a waiver of any or all intellectual property rights that Applied Medical Corporate Group has established in any of its trademarks.

Trademark Name	European Registered
AeroStat	TM
Alexis	®
Applied	TM
Applied Medical	®
Applied Medical Resources	TM
Applied Medical Resources Logo Design	®
Epix	®
Fios	®
Forté	®
GelPOINT	®
GelPort	®
GelSeal	®
Inzii	®
Kii	®
<i>latis</i>	TM
Silhouette	®
Stealth	®
Sureseal	TM
UroWIRE	®
Voyant	®

MADE IN THE USA AND/OR THE NETHERLANDS

© 2023 Applied Medical Resources Corporation. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

Symbol ™ is a reference to a trademark of Applied Medical that has not been registered.

11/2023 196381-EN-N

TABLE OF CONTENTS

ENERGY DEVICES	9
LAPAROSCOPIC DEVICES	
Voyant Laparoscopic Devices	9
OPEN SURGERY DEVICES	
Voyant Open Surgery Devices	9
GENERATOR AND CART	
Voyant Generator and Cart	11
GelSeal TECHNOLOGY	15
HAND ACCESS SYSTEM	
GelPort Laparoscopic System	15
SINGLE-SITE DEVICE	
GelPOINT Advanced Access Platform	17
REDUCED-PORT DEVICE	
GelPOINT Mini Advanced Access Platform	19
TRANSANAL ACCESS PLATFORMS	
GelPOINT Path Transanal Access Platforms	21
TRANSVAGINAL ACCESS PLATFORMS	
GelPOINT V-Path Transvaginal Access Platforms	23
LAPAROSCOPIC SYSTEMS AND DEVICES	27
ABDOMINAL ACCESS SYSTEMS	
Kii Fios First Entry Access Systems	27
Kii Fios First Entry Access System Dual Packs	29
Kii Balloon Blunt-Tip Access Systems	31
Kii Optical Access Systems	33
Kii Optical Access System Dual Packs	35
Kii Shielded Bladed Access Systems	37
Kii Shielded Bladed Access System Dual Packs	39
Kii Sleeves	41
Kii Low-Profile Optical Access Systems	43
Kii Low-Profile Bladed Access Systems	43
Kii Low-Profile Optical Access System Dual Packs	45
Kii Low-Profile Bladed Access System Dual Packs	45
Kii Low-Profile Sleeves	45
Kii Dissecting Balloon Access Systems	47

LAPAROSCOPIC INSTRUMENTATION

Epix Universal Clip Applier	49
Epix Electrosurgical Probes with Smoke Evacuation	51
Epix Laparoscopic Dissector	53
Epix Laparoscopic Graspers	53
Epix Laparoscopic Scissors	53

SUCTION-IRRIGATION SYSTEM

Epix Disposable Tubing Sets	55
Epix Reusable Probes	55

RETRIEVAL SYSTEMS

Inzii Retrieval Systems	57
-------------------------	----

ANCILLARY DEVICES

Applied Scope Warmer	59
Applied Insufflation Needles	59
AeroStat Cholangiography Catheter	61

PROTECTORS AND SYSTEMS 65

GENERAL SURGERY PROTECTORS

Alexis O Wound Protector-Retractors	65
Alexis Wound Protector-Retractors	67
Alexis Laparoscopic Systems	69

CESAREAN SECTION RETRACTORS

Alexis O C-Section Protector-Retractors	71
-----------------------------------------	----

ORTHOPAEDIC PROTECTORS

Alexis Orthopaedic Protectors	73
-------------------------------	----

TISSUE EXTRACTION SYSTEMS

Alexis Contained Extraction Systems	75
Alexis Contained Extraction Systems with Kii Balloon Blunt-Tip Access System	75

ENDOUROLOGY PRODUCTS 79

DEVICES FOR STONE MANAGEMENT

Silhouette Pediatric Ureteral Stents	79
Silhouette Comfort Ureteral Stents	79
Silhouette XtraFlo Ureteral Stents	81
Forté A _x P Ureteral Access Sheaths	83

ENDOUROLOGIC ACCESSORIES

Sureseal II Endoscopic Valve	85
------------------------------	----

VASCULAR AND CARDIOVASCULAR PRODUCTS 89

OCCLUSSION PRODUCTS



Stealth Surgical Clip and Applicator	89
Stealth Surgical Clamp Inserts	91
Stealth Surgical Clamps	93

ENERGY DEVICES

LAPAROSCOPIC DEVICES

OPEN SURGERY DEVICES

GENERATOR AND CART



ENERGY DEVICES

LAPAROSCOPIC DEVICES

Voyant Laparoscopic Devices

Description	Qty/Box	Shaft Dimensions	Reorder Number
Voyant 5mm Fusion device	6	5mm x 37cm	EB210
Voyant 5mm Fusion device	6	5mm x 44cm	EB211
Voyant Maryland Fusion device	6	5mm x 37cm	EB215
Voyant Maryland Fusion device	6	5mm x 44cm	EB216
Voyant Maryland Fusion device	6	5mm x 23cm	EB217

Note: Voyant laparoscopic devices are indicated for both laparoscopic and open procedures.

OPEN SURGERY DEVICES

Voyant Open Surgery Devices

Description	Qty/Box	Shaft Length	Reorder Number
Voyant Fine Fusion device	6	19.3cm ^a	EB230
Voyant Open Fusion device	6	20cm	EB240

^aDevice length.

GENERATOR AND CART



Voyant Generator and Cart

Description	Qty/Box	Reorder Number
Voyant electro-surgical generator	1	EA020
Voyant cart	1	EX150

GelSeal TECHNOLOGY

HAND ACCESS SYSTEM

SINGLE-SITE DEVICE

REDUCED-PORT DEVICE


TRANSANAL ACCESS PLATFORMS

TRANSVAGINAL ACCESS PLATFORMS

GelSeal TECHNOLOGY
HAND ACCESS SYSTEM

GelPort Laparoscopic System

Description	Qty/Box	Reorder Number
GelPort laparoscopic system (1) GelSeal cap (1) Alexis O wound protector-retractor	1	C8XX2CE



GelSeal TECHNOLOGY SINGLE-SITE DEVICE



GelPOINT Advanced Access Platform

Description	Qty/Box	Reorder Number
GelPOINT advanced access platform (1) GelSeal cap (1) Alexis O wound protector-retractor (3) 10mm sleeves (1) 12mm sleeve (1) Obturator for sleeves (1) Instrument shield	1	CNGL2



GelSeal TECHNOLOGY REDUCED-PORT DEVICE

GelPOINT Mini Advanced Access Platform

Description	Qty/Box	Reorder Number
GelPOINT Mini advanced access platform (1) GelSeal cap (1) Alexis O wound protector-retractor (3) 10mm sleeves (1) 12mm sleeve (1) Obturator for sleeves	1	CNGL3



TRANSANAL ACCESS PLATFORMS



4cm access channel
(included in unit CNB10)



5.5cm access channel
(included in unit CNB11)




9cm access channel
(included in unit CNB12)

GelPOINT Path Transanal Access Platforms

Description	Qty/Box	Reorder Number
GelPOINT Path transanal access platform (1) GelSeal cap (1) 4cm access channel with introducer (3) 10mm sleeves (1) 10mm insufflation stabilization bag (ISB) sleeve (1) Obturator for sleeves (1) Insufflation stabilization bag (ISB)	1	CNB10
GelPOINT Path transanal access platform (1) GelSeal cap (1) 5.5cm access channel with introducer (3) 10mm sleeves (1) 10mm insufflation stabilization bag (ISB) sleeve (1) Obturator for sleeves (1) Insufflation stabilization bag (ISB)	1	CNB11
GelPOINT Path transanal access platform (1) GelSeal cap (1) 9cm access channel with introducer (3) 10mm sleeves (1) 10mm insufflation stabilization bag (ISB) sleeve (1) Obturator for sleeves (1) Insufflation stabilization bag (ISB)	1	CNB12

Note: Each GelPOINT Path platform comes with an access channel in one of three available sizes.



GelSeal TECHNOLOGY

TRANSVAGINAL ACCESS PLATFORMS

GelPOINT V-Path Transvaginal Access Platforms

Description	Alexis O Retractor Inner Ring Size	Qty/Box	Reorder Number
GelPOINT V-Path transvaginal access platform (1) GelSeal cap (1) Alexis O retractor (3) 10mm sleeves (1) 12mm sleeve (1) Obturator for sleeves (1) Instrument shield (1) Introducer for the Alexis O retractor	11cm	1	C2A11
GelPOINT V-Path transvaginal access platform (1) GelSeal cap (1) Alexis O retractor (3) 10mm sleeves (1) 12mm sleeve (1) Obturator for sleeves (1) Instrument shield (1) Introducer for the Alexis O retractor	9.5cm	1	C2A12
GelPOINT V-Path transvaginal access platform (1) GelSeal cap (1) Alexis O retractor (3) 10mm sleeves (1) 12mm sleeve (1) Obturator for sleeves (1) Instrument shield (1) Introducer for the Alexis O retractor	7cm	1	C2A15

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS

DISSECTING BALLOON ACCESS SYSTEMS

LAPAROSCOPIC INSTRUMENTATION

SUCTION-IRRIGATION SYSTEM

RETRIEVAL SYSTEMS

ANCILLARY DEVICES

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS



Kii Fios First Entry Access Systems

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 75mm Kii Fios first entry access system	6	CFF05	N/A
5 x 100mm Kii Fios first entry access system	6	CFF03	CTF03
5 x 100mm Kii Fios first entry access system, handled	6	CFF04	CTF04
5 x 150mm Kii Fios first entry access system	6	CFF01	CTF01
11 x 100mm Kii Fios first entry access system	6	CFF33	CTF33
11 x 100mm Kii Fios first entry access system, handled	6	CFF34	CTF34
12 x 75mm Kii Fios first entry access system	6	CFF75	N/A
12 x 100mm Kii Fios first entry access system	6	CFF73	CTF73
12 x 100mm Kii Fios first entry access system, handled	6	CFF74	CTF74
12 x 150mm Kii Fios first entry access system	6	CFF71	CTF71

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS



Kii Fios First Entry Access System Dual Packs

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 75mm Kii Fios first entry access system dual pack (2) Kii sleeves (1) Kii Fios obturator	10	CFF14	N/A
5 x 100mm Kii Fios first entry access system dual pack (2) Kii sleeves (1) Kii Fios obturator	10	CFF12	CTF12
5 x 150mm Kii Fios first entry access system dual pack (2) Kii sleeves (1) Kii Fios obturator	10	CFF10	N/A
11 x 100mm Kii Fios first entry access system dual pack (2) Kii sleeves (1) Kii Fios obturator	10	N/A	CTF35
12 x 100mm Kii Fios first entry access system dual pack (2) Kii sleeves (1) Kii Fios obturator	10	N/A	CTF75

LAPAROSCOPIC SYSTEMS AND DEVICES
ABDOMINAL ACCESS SYSTEMS



Kii Balloon Blunt-Tip Access Systems

Description	Qty/Box	Reorder Number
5 x 100mm Kii balloon blunt-tip access system	6	C0Q20
12 x 100mm Kii balloon blunt-tip access system	6	C0R47
12 x 130mm Kii balloon blunt-tip access system	6	C0R50

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS



Kii Optical Access Systems

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 55mm Kii optical access system	6	N/A	CTR05
8 x 100mm Kii optical access system	6	N/A	COQ19
12 x 60mm Kii optical access system	6	N/A	COR83
15 x 100mm Kii optical access system	6	N/A	COR37
15 x 150mm Kii optical access system	6	CFR39	COR39

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS



Kii Optical Access System Dual Packs

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 55mm Kii optical access system dual pack (2) Kii sleeves (1) Kii optical separator obturator	10	N/A	CTR14
5 x 150mm Kii optical access system dual pack (2) Kii sleeves (1) Kii optical separator obturator	10	CFR10	CTR10



LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS

Kii Shielded Bladed Access Systems

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 55mm Kii shielded bladed access system	6	N/A	CTB05
5 x 100mm Kii shielded bladed access system	6	CFB03	CTB03
5 x 150mm Kii shielded bladed access system	6	CFB01	CTB01
8 x 100mm Kii shielded bladed access system	6	N/A	COQ61
11 x 100mm Kii shielded bladed access system	6	CFB33	CTB33
12 x 100mm Kii shielded bladed access system	6	CFB73	CTB73
12 x 150mm Kii shielded bladed access system	6	CFB71	CTB71

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS



Kii Shielded Bladed Access System Dual Packs

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 55mm Kii shielded bladed access system dual pack (2) Kii sleeves (1) Kii shielded bladed obturator	10	N/A	CTB14
5 x 100mm Kii shielded bladed access system dual pack (2) Kii sleeves (1) Kii shielded bladed obturator	10	CFB12	CTB12
5 x 150mm Kii shielded bladed access system dual pack (2) Kii sleeves (1) Kii shielded bladed obturator	10	N/A	CTB10
11 x 100mm shielded bladed access system dual pack (2) Kii sleeves (1) Kii shielded bladed obturator	10	N/A	CTB35
12 x 100mm Kii shielded bladed access system dual pack (2) Kii sleeves (1) Kii shielded bladed obturator	10	N/A	CTB75

LAPAROSCOPIC SYSTEMS AND DEVICES


ABDOMINAL ACCESS SYSTEMS



Kii Sleeves

Description	Qty/Box	Advanced Fixation Reorder Number	Z-Thread Reorder Number
5 x 55mm Kii sleeve	12	N/A	CTS03
5 x 75mm Kii sleeve ^a	12	CFS03	N/A
5 x 100mm Kii sleeve	12	CFS02	CTS02
5 x 150mm Kii sleeve	12	CFS01	CTS01
11 x 100mm Kii sleeve	12	CFS12	CTS12
12 x 100mm Kii sleeve	12	CFS22	CTS22
12 x 150mm Kii sleeve	12	CFS21	CTS21

^a For use with the Kii Fios obturator, 5 x 75mm.



LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS

Kii Low-Profile Optical Access Systems


Description	Qty/Box	Reorder Number
5 x 55mm Kii low-profile optical access system	6	CTR21
5 x 100mm Kii low-profile optical access system	6	CTR22

Kii Low-Profile Bladed Access Systems

Description	Qty/Box	Reorder Number
5 x 55mm Kii low-profile bladed access system	6	CTB21
5 x 100mm Kii low-profile bladed access system	6	CTB22

LAPAROSCOPIC SYSTEMS AND DEVICES

ABDOMINAL ACCESS SYSTEMS



Kii Low-Profile Optical Access System Dual Packs

Description	Qty/Box	Reorder Number
5 x 55mm Kii low-profile optical access system dual pack (2) Kii low-profile sleeves (1) Kii low-profile optical separator obturator	10	CTR23
5 x 100mm Kii low-profile optical access system dual pack (2) Kii low-profile sleeves (1) Kii low-profile optical separator obturator	10	CTR24

Kii Low-Profile Bladed Access System Dual Packs

Description	Qty/Box	Reorder Number
5 x 55mm Kii low-profile bladed access system dual pack (2) Kii low-profile sleeves (1) Kii low-profile bladed obturator	10	CTB23
5 x 100mm Kii low-profile bladed dual pack (2) Kii low-profile sleeves (1) Kii low-profile bladed obturator	10	CTB24

Kii Low-Profile Sleeves

Description	Qty/Box	Reorder Number
5 x 55mm Kii low-profile sleeve ^a	12	CTS28
5 x 100mm Kii low-profile sleeve ^a	12	CTS29

^aCompatible only with Kii low-profile obturators.

LAPAROSCOPIC SYSTEMS AND DEVICES
**DISSECTING BALLOON
ACCESS SYSTEMS**



Kii Dissecting Balloon Access Systems

Description	Qty/Box	Reorder No.
Kii dissecting balloon access system, oval	3	C0H11
Kii dissecting balloon access system, round	3	C0H12



LAPAROSCOPIC SYSTEMS AND DEVICES
LAPAROSCOPIC INSTRUMENTATION

Epix Universal Clip Applier

Description	Qty/Box	Clip Size	Port Size	Reorder Number
Epix universal clip applier	3	M/L	5mm ^a	CA500

^aApplied Medical Kii trocar (if using non-Kii trocars, a larger port size will be needed).



LAPAROSCOPIC SYSTEMS AND DEVICES

LAPAROSCOPIC INSTRUMENTATION



Epix Electrosurgical Probes with Smoke Evacuation

Description	Qty/Box	Overall Length	Reorder Number
Probe with angled L-hook tip	1	5mm x 42cm ^a	CW002
Probe with straight L-hook tip	1	5mm x 42cm ^b	CW005

^a Compatible with Applied Medical GelPOINT sleeves and Applied Medical Kii 12 x 60-100mm and 15 x 100-150mm trocars.

^b Compatible with Applied Medical GelPOINT sleeves and Applied Medical Kii 5mm trocars (if using non-Kii trocars, a larger port size will be needed).



LAPAROSCOPIC SYSTEMS AND DEVICES

LAPAROSCOPIC INSTRUMENTATION



Epix Laparoscopic Dissector

Description	Qty/Box	Overall Length	Reorder Number
Epix laparoscopic dissector	10	5 x 35cm	CY010

Epix Laparoscopic Graspers

Description	Qty/Box	Overall Length	Reorder Number
Epix <i>latis</i> graspers	10	5 x 35cm	C4130
	10	5 x 45cm	C4140

Epix Laparoscopic Scissors

Description	Qty/Box	Overall Length	Reorder Number
Epix laparoscopic scissors	10	5 x 35cm	CB030
	10	5 x 45cm	CB040

LAPAROSCOPIC SYSTEMS AND DEVICES

SUCTION-IRRIGATION SYSTEM



Epix Disposable Tubing Sets

Description	Qty/Box	Diameter	Working Length	Reorder Number
Probe and valve with dual-bag tubing	12	5mm	33cm	C6001
Valve only with dual-bag tubing	12	N/A	N/A	C6003

Epix Reusable Probes

Description	Qty/Box	Diameter	Working Length	Reorder Number
Probes with irrigation holes	1	5mm	45cm	C7211
	1	10mm	33cm	C7204



LAPAROSCOPIC SYSTEMS AND DEVICES

RETRIEVAL SYSTEMS

Inzii Retrieval Systems

Description	Qty/Box	Port Size	Volume	Reorder Number
Inzii universal retrieval system	10	5mm ^a	180mL	CD003
Inzii 10mm retrieval system	10	10mm	225mL	CD001
Inzii 12/15mm retrieval system	5	12mm	1600mL	CD004

^aApplied Medical Kii trocar (if using a non-Kii trocar, a larger port size will be needed).



LAPAROSCOPIC SYSTEMS AND DEVICES

ANCILLARY DEVICES



Applied Scope Warmer

Description	Qty/Box	Reorder Number
Disposable scope-warmer seal	20	C3101
Reusable scope-warmer canister ^a	1	C3001
Reusable scope-warmer base ^a	1	C3002

^aNonsterile.


Applied Insufflation Needles

Description	Qty/Box	Reorder Number
Applied insufflation needle, 13 gauge, 120mm, disposable	20	C2201
Applied insufflation needle, 13 gauge, 150mm, disposable	20	C2202

LAPAROSCOPIC SYSTEMS AND DEVICES
ANCILLARY DEVICES

AeroStat Cholangiography Catheter

Description	Qty/Box	Reorder Number
AeroStat cholangiography catheter, 3F tip, with percutaneous introducer, 6-20F retention cone, disposable	5	C1002



PROTECTORS AND SYSTEMS

GENERAL SURGERY PROTECTORS

CESAREAN SECTION RETRACTORS

ORTHOPAEDIC PROTECTORS

TISSUE EXTRACTION SYSTEMS

GENERAL SURGERY PROTECTORS



Alexis O Wound Protector-Retractors

Featuring a rigid retraction ring

Description	Qty/Box	Incision Range	Sheath Length	Reorder Number
Small ^a	5	2.5–6cm	18cm	C8401
Medium	5	5–9cm	18cm	C8402
Large	5	9–14cm	25cm	C8403
Extra large	5	11–17cm	34cm	C8404
Extra extra large	5	17–25cm	36cm	C8405
Extra extra extra large	3	25–32cm	39cm	C8406

^aModel includes a tether to facilitate device removal.

GENERAL SURGERY PROTECTORS



Alexis Wound Protector-Retractors

Featuring a flexible retraction ring

Description	Qty/Box	Incision Range	Sheath Length	Reorder Number
Extra extra small ^a	5	1–3cm	20cm	C8313
Extra extra small, short ^a	5	1–3cm	11cm	C8323
Extra small ^a	5	2–4cm	19cm	C8312
Extra small, short ^a	5	2–4cm	13cm	C8322
Small ^a	5	2.5–6cm	18cm	C8301
Medium	5	5–9cm	18cm	C8302
Large	5	9–14cm	25cm	C8303
Extra large	5	11–17cm	34cm	C8304

^a Model includes a tether to facilitate device removal.

GENERAL SURGERY PROTECTORS



Alexis Laparoscopic Systems

Featuring a laparoscopic cap

Description	Qty/Box	Incision Range	Sheath Length	Reorder Number
Alexis laparoscopic system (1) Alexis laparoscopic cap, small (1) Alexis O wound protector-retractor, small	6	2.5–6cm	18cm	C8501 ^a
Alexis laparoscopic system (1) Alexis laparoscopic cap, medium (1) Alexis O wound protector-retractor, medium	6	5–9cm	18cm	C8502

^a Model includes a tether to facilitate device removal.

CESAREAN SECTION RETRACTORS



Alexis O C-Section Protector-Retractors

Featuring a rigid retraction ring

Description	Qty/Box	Incision Range	Sheath Length	Reorder Number
Large	5	9-14cm	25cm	G6313
Extra large	5	11-17cm	34cm	G6314

ORTHOPAEDIC PROTECTORS



Alexis Orthopaedic Protectors

Featuring a rigid retraction ring

Description	Qty/Box	Incision Range	Sheath Length	Reorder Number
Extra small/medium	5	2.5–7cm	14cm	HR000
Small/small	5	2.5–8cm	14cm	HR001
Small/medium	5	2.5–8cm	14cm	HR004
Medium/large	5	5–13cm	17cm	HR005

Alexis Orthopaedic Protectors

Featuring a flexible retraction ring

Description	Qty/Box	Incision Range	Sheath Length	Reorder Number
Extra small/medium	5	2.5–7cm	14cm	HR100
Small/small	5	2.5–8cm	14cm	HR101
Small/medium	5	2.5–8cm	14cm	HR104
Medium/large	5	5–13cm	17cm	HR105

TISSUE EXTRACTION SYSTEMS



Alexis Contained Extraction Systems

Description	Qty/ Box	Ring Diameter	Volume	Reorder Number
Alexis contained extraction system, 14cm (1) Specimen containment bag, 14cm diameter (1) Small guard (1) Medium guard (1) Retractor	3	14cm	3400mL	GTB14
Alexis contained extraction system, 17cm (1) Specimen containment bag, 17cm diameter (1) Small guard (1) Medium guard (1) Retractor	3	17cm	6500mL	GTB17


Alexis Contained Extraction Systems with Kii Balloon Blunt-Tip Access System

Description	Qty/ Box	Ring Diameter	Volume	Reorder Number
Alexis contained extraction system, 14cm, with Kii balloon blunt-tip access system (1) Specimen containment bag, 14cm diameter (1) Small guard (1) Medium guard (1) Retractor (1) 12 x 100mm Kii balloon blunt-tip access system	1	14cm	3400mL	GK1059
Alexis contained extraction system, 17cm, with Kii balloon blunt-tip access system (1) Specimen containment bag, 17cm diameter (1) Small guard (1) Medium guard (1) Retractor (1) 12 x 100mm Kii balloon blunt-tip access system	1	17cm	6500mL	GK1060

Please contact your field representative for kit availability.

ENDOUROLOGY PRODUCTS

DEVICES FOR STONE MANAGEMENT
ENDOUROLOGIC ACCESSORIES



ENDOUROLOGY PRODUCTS

DEVICES FOR STONE MANAGEMENT

Silhouette Pediatric Ureteral Stents

Description	Qty/Box	Size	Length	Reorder Number
Silhouette Pediatric ureteral stents with 0.035" UroWIRE guidewire	1	4F	10cm	B3884

Silhouette Comfort Ureteral Stents

Description	Qty/Box	Size	Length	Reorder Number
Silhouette Comfort ureteral stents	1	4.6F	24cm	B3836
	1	4.6F	26cm	B3837
	1	4.6F	28cm	B3838

ENDOUROLOGY PRODUCTS
DEVICES FOR STONE MANAGEMENT

Silhouette XtraFlo Ureteral Stents

Description	Qty/Box	Size	Length	Reorder Number
Silhouette XtraFlo ureteral stents	1	6F	24cm	B3856
	1	6F	26cm	B3857
	1	6F	28cm	B3858
Silhouette XtraFlo multi-length stent	1	6F	22–30cm	B3847



DEVICES FOR STONE MANAGEMENT



Forté A_xP Ureteral Access Sheaths

Description	Qty/Box	Size		Length	Reorder Number
		OD Dilator	OD Sheath		
Size 1	6	10F	12F	28cm	B7276
	6	10F	12F	35cm	B7216
	6	10F	12F	45cm	B7296
	6	10F	12F	55cm	B7286
Size 2	6	12F	14F	20cm	B7066
	6	12F	14F	35cm	B7016
	6	12F	14F	45cm	B7096
	6	12F	14F	55cm	B7086
Size 3	6	14F	16F	20cm	B7166
	6	14F	16F	28cm	B7176
	6	14F	16F	35cm	B7116

ENDOUROLOGIC ACCESSORIES



Sureseal II Endoscopic Valve

Description	Qty/Box	Size	Reorder Number
Sureseal II endoscopic valve	10	0-5F	SEAL-005

VASCULAR AND CARDIOVASCULAR PRODUCTS

OCCLUSION PRODUCTS

VASCULAR AND CARDIOVASCULAR PRODUCTS

OCCLUSION PRODUCTS



Stealth Surgical Clip and Applicator

For the IMA, saphenous vein, renal artery and deep femoral artery

Description	Qty/Box	Size	Color	Reorder Number
Stealth <i>latis/latis</i> clips	10	6mm (1/2 force)	Dark blue	A1601
	10	6mm (1/4 force)	Light blue	A1602
	10	6mm (3/4 force)	White	A1603
Stealth Soft/Soft clip	10	6mm (3/4 force)	White	A1604
Stealth surgical clip applicator	1			A3200

latis/latis clips



Soft/Soft clip



OCCLUSION PRODUCTS





Stealth Surgical Clamp Inserts

Description	Qty/Box	Size	Reorder Number
Stealth <i>latis/latis</i>	10	30mm (#1)	A0C01
	10	45mm (#2)	A0C02
	10	60mm (#3)	A0C03
	10	90mm (#5)	A0C04
	10	Small sidebite	A0C05
	10	Medium sidebite	A0C06

latis/latis inserts



OCCLUSION PRODUCTS



Stealth Surgical Clamps

Applications	Handle Angle	Overall Length	Jaw Angle	Jaw Length	Jaw Size	Reorder Number
Femoral and iliac	15°	16.5cm	Straight	30mm	Size 1	A3201
	25°	15.2cm	60°	30mm	Size 1	A3202
	15°	17.8cm	Straight	45mm	Size 2	A3203
	15°	17.8cm	34°	45mm	Size 2	A3204
	Straight	16.5cm	65°	45mm	Size 2	A3205
	15°	16.5cm	65°	45mm	Size 2	A3206
Iliac and abdominal aortic	10°	29.2cm	Straight	60mm	Size 3	A3207
	15°	25.4cm	60°	60mm	Size 3	A3208
	35°	25.4cm	60°	60mm	Size 3	A3209
	35°	25.4cm	Bahnson	60mm	Size 3	A3210
	Iliac, abdominal aortic, femoral and renal	60°	26.7cm	Wylie	60mm	Size 3
Aortic cross-clamp	15°	30.5cm	Straight	90mm	Size 5	A3212
	35°	27.9cm	60°	90mm	Size 5	A3213
	35°	26.7cm	Bahnson	90mm	Size 5	A3214
	35°	22.9cm	Bahnson	90mm	Size 5	A3220
Aortic and vena cava	25°	17.8cm	Sidebite	50mm	Small	A3215
	45°	20.3cm	Sidebite	75mm	Medium	A3216

