




Stealth[®]

Atraumatic Occlusion System

Advantages of Stealth Clips

Stealth clips are designed to provide occlusion where space is at a premium.



Low Profile

Ultrasoft Pad



Reduced Spring Force

Three options allow for 1/4-force, 1/2-force or 3/4-force occlusion, depending on the surgeon's preference and the needs of each procedure.



latis[®]/*latis* Clips

Reorder No.	Description	Size
A1601	<i>latis/latis</i> clip	6mm (1/2 Force)
A1602	<i>latis/latis</i> clip	6mm (1/4 Force)
A1603	<i>latis/latis</i> clip	6mm (3/4 Force)



Soft/Soft Clips

Reorder No.	Description	Size
A1604	Soft/Soft clip	6mm (3/4 Force)



Stealth Surgical Clip Applicator

Reorder No.	Description	Handle Angle	Overall Length
A3200	Stealth surgical clip applicator	13°	18.4cm

latis Inserts

The unique *latis* mesh insert gently engages the adventitia reducing migration and compressive force.




latis/latis Inserts

Reorder No.	Size
<i>Straight</i>	
A0C01	30mm
A0C02	45mm
A0C03	60mm
A0C04	90mm
<i>Sidebites</i>	
A0C05	Small
A0C06	Medium

Advantages of the Stealth Occlusion System

The low-profile Stealth system has become a leader in atraumatic occlusion because of its excellent traction, reduced occlusion force and minimized clip migration.

Stealth Clamps Versus Standard Clamps



Ease of Installation

The inserts attach securely and easily.



Variety of Sizes

Stealth systems offer a wide variety of clamp insert sizes and configurations.

Size 1 - 30mm




Size 2 - 45mm




Size 3 - 60mm



Size 5 - 90mm



Sidebite, small
50mm



Sidebite, medium
75mm

Actual sizes

Stealth Clamps

Atraumatic Stealth clamps offer an unmatched level of versatility.



Ensured Durability


Stainless steel construction ensures durability.

Enhanced Stability

The exceptionally crafted, fine-tuned box lock enhances stability.

Consistent, Gentle Occlusion

The jaws provide consistent, gentle occlusion as a result of their parallel alignment.



Low Profile

The low-profile design enables optimal occlusion in hard-to-access areas.

Firmly Secured Inserts

Innovative attachment technology fastens the insert securely to the clamp.

Preferred Configurations



Bahnsen

The Stealth system's best-selling model



Wylie

Deep reach with a unique bend



Sidebites

The market's only padded sidebite

Stealth Clamp Configurations



Reorder no.	Size	Handle Angle	Overall Length	Jaw Angle	Jaw Length	Applications
A3201	1	15°	16.5cm	Straight	30mm	Femoral and iliac
A3202	1	25°	15.2cm	60°	30mm	Femoral and iliac
A3203	2	15°	17.8cm	Straight	45mm	Femoral and iliac
A3204	2	15°	17.8cm	34°	45mm	Femoral and iliac
A3205	2	Straight	16.5cm	65°	45mm	Femoral and iliac
A3206	2	15°	16.5cm	65°	45mm	Femoral and iliac
A3207	3	10°	29.2cm	Straight	60mm	Iliac and abdominal aortic
A3208	3	15°	25.4cm	60°	60mm	Iliac and abdominal aortic

Stealth Clamp Configurations



Reorder no.	Size	Handle Angle	Overall Length	Jaw Angle	Jaw Length	Applications
A3209	3	35°	25.4cm	60°	60mm	Iliac and abdominal aortic
A3210	3	35°	25.4cm	Bahnsen	60mm	Iliac and abdominal aortic
A3211	3	60°	26.7cm	Wylie	60mm	Iliac, abdominal aortic, femoral and renal
A3212	5	15°	30.5cm	Straight	90mm	Aortic cross-clamp
A3213	5	35°	27.9cm	60°	90mm	Aortic cross-clamp
A3214	5	35°	26.7cm	Bahnsen	90mm	Aortic cross-clamp
A3215	Small	25°	17.8cm	Sidebite	50mm	Aortic and vena cava
A3216	Medium	45°	20.3cm	Sidebite	75mm	Aortic and vena cava
A3220	5	35°	22.9cm	Bahnsen	90mm	Aortic cross-clamp

Visit [appliedmedical.com/stealth](https://www.appliedmedical.com/stealth) for more information.

Stealth clips, clamps and inserts are used for surgical clamping during cardiovascular, peripheral vascular and general surgery.

Please contact your Applied Medical representative for more information on availability. This information is intended for dissemination exclusively to healthcare professionals and is not intended to replace labeling and Instructions for Use (IFU). Please refer to the IFU for the indications, contraindications, warnings, precautions, instructions and other information.

© 2022 Applied Medical Resources Corporation, Rancho Santa Margarita, CA. All rights reserved. Applied Medical, the Applied Medical logo design, and marks designated with the symbol ® are trademarks of Applied Medical Resources Corporation, registered in one or more of the following countries: Australia, Canada, Japan, South Korea, the United States, the United Kingdom, and/or the European Union.

08/2022 271883-EN-B

