


Sureseal[®] II

Endoscopic Valve



The background of the slide features a large, abstract teal shape on the left side, composed of several overlapping triangles and polygons, creating a modern, geometric look. The rest of the background is white.

THE SEAL THAT KEEPS YOU DRY

The Sureseal II endoscopic valve combines maximum irrigation capability with a leak-tight seal around 0–5F stone baskets, guidewires, laser fibers and other endourological instruments.

Superior Closure

Keeps urologist dry

Works with a variety of scopes and the Applied Forté ureteral access sheaths

Accommodates instruments from 0-5F

Enhanced Handling

Allows passage of instruments into working port without additional manipulation of valve

Enables one-handed manipulation of scope

Latex Free

Protects latex sensitive patients





Sureseal II Endoscopic Valve

Reorder No.	Valve Opening	Pressure Rating	Irrigation Port	Scope Port	Qty/Box
Seal-005	Passive, continuously reliable bidirectional from 0-5F (.065")	20psi or 1034mm Hg	Luer-lock type fitting	Standard luer-lock with ball adapter	10/BOX (includes valve and ball adapter)

Visit www.appliedmedical.com/accessories
for more information